



INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification ⁶ : G01N 33/00, 33/48	A1	(11) International Publication Number: WO 99/58968 (43) International Publication Date: 18 November 1999 (18.11.99)
(21) International Application Number: PCT/IL99/00252 (22) International Filing Date: 12 May 1999 (12.05.99) (30) Priority Data: 124444 12 May 1998 (12.05.98) IL (71) Applicant (for all designated States except US): ADVANCED NEUROPROTECTIVE SYSTEMS LTD. [IL/IL]; Kiryat Weizman Science Park, 70400 Nes Ziona (IL). (72) Inventors; and (75) Inventors/Applicants (for US only): DOLINA, Svetlana [IL/IL]; Aharoni 14/12, 76282 Rehovot (IL). RABINKOV, Aaron [IL/IL]; Bnei Moshe 19/8, 76485 Rehovot (IL). PRESSMAN, Eugene [IL/IL]; Hadadi Zion Shaul 21/22, 59514 Bat Yam (IL). (74) Agent: FRIEDMAN, Mark, M.; Beit Samueloff Building, Haomanim 7, 67897 Tel Aviv (IL).		(81) Designated States: AE, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CU, CZ, DE, DK, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, UA, UG, US, UZ, VN, YU, ZA, ZW, ARIPO patent (GH, GM, KE, LS, MW, SD, SL, SZ, UG, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG). Published <i>With international search report.</i> <i>With amended claims.</i>
(54) Title: DIAGNOSIS OF PREDISPOSITION TO EPILEPSY AND MONITORING OF ANTIEPILEPTIC TREATMENT (57) Abstract <p>A method and a diagnostic system for the detection of a predisposition to epilepsy in a patient, for the diagnosis of clinical epilepsy itself and for monitoring of antiepileptic treatment. The method includes the steps of obtaining a sample in a patient, preferably a blood sample, determining a concentration in the sample of at least two kynurenine metabolites selected from the group consisting of TRP, L-KYN, KA, 3HOAA, AA and QUIN, those having functionally opposite properties, and correlating the concentration of these metabolites in the sample and the ratio between them with a range of values of these concentrations and ratios between them in normal individuals. In addition, the level of activity of enzymes in peripheral tissues such as blood cells can also be used to diagnose seizure disorders in the subject, as well as to detect a predisposition to such seizure disorders in the subject. For monitoring of antiepileptic treatment, the range of values of two of the metabolites and the ratio between them are compared with the ranges of values of concentrations and the ratio between them in epileptic patients well controlled or non-controlled by antiepileptic drugs.</p>		